

Alzheimer's Disease



Aβ Seeding

The pathological aggregation of amyloid- β (A β) peptides is one of the major causes for the progressive cognitive decline in Alzheimer's disease patients. The development of compounds that interfere with A β seeding or aggregation and that are thus able to rescue neurodegeneration is indispensable. For this purpose, fast and reliable assays are needed that can show direct effects of compounds on A β oligomer formation, seeding and aggregation.

To screen for developmental compounds, interfering with A β seeding and aggregation properties, kinetics of A β aggregation can be monitored over time by using Thioflavin T dye (ThT). The assay is based on the property of ThT dye in which fluorescence (Ex/Em=440/484 nm) is increased when bound to aggregated A β peptides. In addition, QPS Neuropharmacology exclusively provides a fast and reproducible screening assay, Amorfix Aggregated A β Assay (A4), that can specifically assess beneficial effects of your developmental compounds on the formation of A β oligomers.

Figure 1: Determination of $\ensuremath{\mathsf{A}\beta}$ aggregation and seeding in vitro by ThT Assay. Compounds were co-aggregated with A_β1-42 and ThT and kinetics of Aß aggregation are monitored over time as increase in green fluorescence. Data are presented as RFU (relative fluorescence unit). The reference item at both concentrations (A; RI, tannic acid) as well as the test item (B) at 100 µM, but not 12.5 uM. concentration were found to reduce Aβ aggregation. Mean ± SD;

Figure 2: Determination of Aß oligomers in vitro by A4 Assay. Compounds were co-aggregated with Aβ1-42 for 48 h in vitro. Aggregated Aß was separated from monomers through affinity interaction. After disaggregation, the originally aggregated Aβ was detected using an immunosorbent assay. <u>Aβ</u> levels were evaluated as % aggregated Aß of vehicle control (VC). Mean ± SEM; n=6; One-way ANOVA with Bonferroni's post hoc test; *p<0.001.





Figure 2
A4 assay



Scantox

Discovery

Important note

Representative data are shown throughout this document. However, biological variability might cause deviations from shown data.

Scantox Group, HQ

Hestehavevej 36A, Ejby
DK – 4623 Lille Skensved
clientservice@scantox.com
www.scantox.com
+45 5686 1500

© Scantox A/S

Scantox is a registered trademark of Scantox A/S.

