



# $\alpha$ -synuclein Pre-Formed Fibril-Induced Aggregation in SNCA and GBA Mutant iPSC-Derived Neurons and Rescue by Antibody Therapy

Loeffler T.<sup>1</sup>, Oosterveen T.<sup>2</sup>, Turner A.<sup>2</sup>, Perez-Ruiz I.<sup>1</sup>, Schilcher I.<sup>1</sup>, Daurer M.<sup>1</sup>, Breznik L.<sup>1</sup>, Flunkert S.<sup>1</sup>, Prokesch M.<sup>1</sup>

<sup>1</sup>Scantox Neuro GmbH, Grambach, Austria | <sup>2</sup>bit.bio Ltd., Cambridge, United Kingdom

**bit.bio**  
THE CELL CODING COMPANY

**scantox**




For further information and inquiries contact office-austria@scantox.com

## Objectives


To investigate how different Parkinson's disease-associated genotypes influence  $\alpha$ -synuclein ( $\alpha$ -syn) pathology and neuronal damage, we compared induced pluripotent stem cell (iPSC)-derived glutamatergic neurons carrying mutations in SNCA (A53T/A53T) and GBA (null/null, null/N409S, null/R159W) with their isogenic controls. We aimed to assess glucocerebrosidase (GCase) activity, neurofilament light chain (NF-L) secretion, and  $\alpha$ -syn levels at baseline and following exposure to  $\alpha$ -syn pre-formed fibrils (PFFs) and therapeutic antibody treatment.

## Results


Neurons harboring GBA mutations displayed reduced GCase activity and elevated NF-L secretion compared to isogenic controls. Baseline  $\alpha$ -syn levels were only slightly increased in mutant lines. Upon PFF treatment, a marked accumulation of both soluble and insoluble  $\alpha$ -syn was observed in all genotypes with higher insoluble level in GBA mutated lines. Antibody treatment effectively reduced both soluble and insoluble  $\alpha$ -syn species across genotypes.



**Figure 1: GCase activity ioGlutamatergic Neurons.** GCase activity was measured on DIV9 on different types of ioGlutamatergic Neurons (25,000 cells/well, 96wp) with SNCA A53T mutations: SNCA A53T/A53T and SNCA A53T/WT as well as GBA gene defects (GBA null/null, GBA null/R159W and GBA null/N409S), wild type (WT) control neurons were also included in the experiment. One-way ANOVA followed by Dunnett's multiple comparison test. \*\*\*p < 0.001.



**Figure 2: NF-L levels in supernatant of ioGlutamatergic Neurons.** NF-L baseline levels were measured at DIV10 in supernatants of different types of ioGlutamatergic Neurons (25,000 cells/well, 96wp) with SNCA A53T mutations: SNCA A53T/A53T as well as GBA gene defects (GBA null/null and GBA null/N409S), wild type (WT) control neurons were also included in the experiment. NF-light® ELISA 10-7001 CE from UmanDiagnostics was used. One-way ANOVA followed by Dunnett's multiple comparison test. \*\*p < 0.01; \*\*\*p < 0.001.



**Figure 3: Soluble (SOL) and insoluble (INSOL) human  $\alpha$ -syn levels in the cell extract from ioGlutamatergic Neurons of different genotypes.** Cells were cultured for 12 days and Triton X-100 soluble and insoluble human  $\alpha$ -syn levels were assessed by MSD assay. Baselines measurement (A, B) and response to  $\alpha$ -syn PFF (C, D) are shown. Data are given as pg  $\alpha$ -syn per  $\mu$ g total protein. Mean + SEM; n = 5-6; One-Way ANOVA followed by Dunnett's multiple comparison test. \*p < 0.05; \*\*p < 0.01; \*\*\*p < 0.001.

## Conclusion


The increased susceptibility of mutant neurons to  $\alpha$ -syn aggregation following PFF challenge underscores their disease relevance. Importantly,  $\alpha$ -syn antibody treatment mitigated fibril-induced  $\alpha$ -syn pathology, supporting therapeutic strategies targeting  $\alpha$ -syn in genetically defined subtypes of Parkinson's disease.

## Materials and Methods

Human iPSC-derived ioGlutamatergic Neurons from mutant and isogenic control lines were cultivated and deterministically programmed according to established protocols by bit.bio. Baseline GCase activity, NF-L release, and  $\alpha$ -syn expression were quantified across genotypes. Cultures were then treated with recombinant  $\alpha$ -syn PFFs to induce aggregation. Soluble and insoluble  $\alpha$ -syn fractions were analyzed via an immunosorbent assay. To test therapeutic modulation, fibril-challenged neurons were treated with an  $\alpha$ -syn antibody, and effects on  $\alpha$ -syn species were quantified with the same methods.


## Results

### Induction of $\alpha$ -syn Aggregation with $\alpha$ -syn PFF



**Figure 4: Soluble and insoluble human  $\alpha$ -syn levels in the cell extract from ioGlutamatergic Neurons of different genotypes.** Wild type (WT) control (A, F), GBA null/null (B, G), GBA null/N409S (C, H), GBA null/R159W (D, I) & SNCA A53T/A53T (E, J) neurons were treated with 0.4  $\mu$ M  $\alpha$ -syn monomers, oligomers or pre-formed fibrils for 12 days and Triton X-100 soluble (SOL) and insoluble (INSOL) human  $\alpha$ -syn levels were assessed by MSD assay. Data are given as pg  $\alpha$ -syn per  $\mu$ g total protein. Mean + SEM; n = 5; One-Way ANOVA followed by Dunnett's multiple comparison test. \*\*p < 0.01; \*\*\*p < 0.001.

### Modulation of $\alpha$ -syn Aggregation with $\alpha$ -syn Antibody



**Figure 5: Soluble and insoluble human  $\alpha$ -syn levels in the cell extract from treated ioGlutamatergic Neurons of different genotypes.** Wild type (WT) control (A, F), GBA null/null (B, G), GBA null/N409S (C, H), GBA null/R159W (D, I) & SNCA A53T/A53T (E, J) neurons were treated with  $\alpha$ -syn pre-formed fibrils ( $\alpha$ -syn PFF) at 0.3  $\mu$ M alone or in combination with anti- $\alpha$ -syn antibody [4D6] at 20 nM for 12 days and Triton X-100 soluble (SOL) and insoluble (INSOL) human  $\alpha$ -syn levels were assessed by MSD assay. Data are given as pg  $\alpha$ -syn per  $\mu$ g total protein. Mean + SEM; n = 5-6; One-Way ANOVA followed by Dunnett's multiple comparison test. \*\*\*p < 0.001.

**Meet Scantox at booth #67**

For more information about the models please visit: [www.scantox.com](http://www.scantox.com)

Scantox Neuro GmbH, Parkring 12, 8074 Grambach, Austria